

Alcohol free Disinfection wipes

For Surfaces and Medical Devices



Effective against NOROVIRUS in  Cleans and disinfects in one

PRODUCT SPECIFICATIONS

Sani-Cloth AF Universal

are disposable disinfection wipes for non-porous hard surfaces and non-invasive medical devices. Sani-Cloth AF Universal wipes can be used to pre-clean or decontaminate medical devices prior to sterilisation or high level disinfection.

Test Method			Contact Time
Bacteria	EN16615	<i>Pseudomonas aeruginosa</i>	60 seconds
	EN16615	<i>Staphylococcus aureus</i>	60 seconds
	EN16615	<i>Enterococcus hirae</i>	60 seconds
	EN13727	<i>Pseudomonas aeruginosa</i>	30 seconds
	EN13727	<i>Staphylococcus aureus</i>	30 seconds
	EN13727	<i>Enterococcus hirae</i>	30 seconds
Enveloped Virus	EN14476	HBV, HCV, HIV, HSNT, SARS, Corona	30 seconds
	EN14476	Vaccinia Virus – Strain Elstree	30 seconds
Non-enveloped Virus	EN14476	Norovirus	30 seconds
Yeasts	EN16615	<i>Candida albicans</i>	60 seconds
	EN13624	<i>Candida albicans</i>	30 seconds
TB	EN14348	<i>Mycobacterium terrae</i>	2 minutes

Multipurpose
Cleans and
disinfects
in one

Broad
spectrum of
kill with low
contact
times

Cost effective
via product
standardisation

EN 16615
Supported by
extracted liquid
testing

Alcohol Free
PHMB Free

Sani-Care
Training
Package*

*Terms & Conditions
Apply

CE
2797

Class IIa
Medical
Device

200 Count Flow Wrap Product code: XP00390 NHS code: VJT592
Wipe size: 267mm x 210mm

40 Count Flow Wrap Product code: XP00391 NHS code: VJT640
Wipe size: 200mm x 200mm

100 Count Flow Wrap Product code: XP00363 NHS code: VJT521
Wipe size: 280mm x 200mm

200 Count Canister Product code: XP00366 NHS code: VJT514
Wipe size: 200mm x 200mm

225 Count Bucket Product code: XP00365 NHS code: VJT513
Wipe size: 270mm x 267mm